
# OBUSFORME®

## COMFORT-PULL™ BACK BRACE

BB-CPB-SM /BB-CPB-MD /BB-CPB-LG

### CONGRATULATIONS

on your purchase of the ObusForme® Back Brace. When used together with a complete back safety training program, your back belt is designed to support your lower back and abdominal muscles to help promote proper posture and good body mechanics while lifting, pushing, pulling, walking, and repetitive bending or twisting.

### INSTRUCTIONS

You can wear your back brace either on top of or under your outer clothing.


1. Separate the two ends of the back brace (at the ObusForme logo) and spread the brace out on a flat surface.
2. Detach the pull handles from the main belt and reattach approximately one or two inches (2.5 or 5 cm) away from the pulley mechanism.

The cords attached to the pull handle will be very loose.

3. Pull the two halves of the back brace belt as far apart as possible so that the cords between the two halves are taut. This will create space to tighten the brace with the Comfort-Pull system after you have fastened the back brace in place.


4. Wrap the back brace around your waist, ensuring that the space between the two halves is centered over your spine.

Overlap the end with the ObusForme logo and handle over the opposite end until the brace feels snug and press the ends together to fasten in place.

5. Using both hands, detach each pull handle from the belt and simultaneously pull them the same distance away from your body until the desired level of compression has been reached. Keeping the cords taut, press the handles against the belt to fasten them into place.

6. To tighten or loosen the back brace, detach the two pull handles and either slacken the pull cords or pull them tighter, re-fastening when the desired compression level has been reached.
7. To ensure optimal results every time you use the back brace, loosen the pull cords and widen the space between the two halves before every use.

Caution: Do not machine wash or dry. Do not bleach. To clean your back brace, spot clean with a wash cloth and mild soapy water; ensure that all soap residue is then removed. For a more comprehensive cleaning, the two pulley mechanisms, with their pull handles, can be detached from the fabric. The two fabric halves of the fabric belt can then be hand washed in cold water with mild detergent. Rinse with cold water to remove all soap residue. Dry with a towel, then lay flat to air dry.

Ensure that the two pulley mechanisms are stored in a safe place during the cleaning process and that the cords between them are not twisted or damaged. Do not submerge the pulley mechanisms in water.

Please Note: Do not permit your back brace to become a substitute for good lifting, bending, pulling, pushing or twisting techniques.

Your back brace will not enable you to lift heavier objects than you normally can. Do not subject yourself to any risk by attempting to lift more weight than you can easily handle.

Only wear your back brace if you are performing tasks or participating in activities where proper posture and good body mechanics are an issue. Continuous wear may cause loss of muscle tone.

### OBUSFORME LIMITED WARRANTY

ObusForme guarantees all items are free from defects in workmanship & materials for a period of time between the original purchase date & that stated below. This guarantee applies when items are used for the purpose intended. Items will be repaired/replaced (at our option), with new, refurbished parts/products &/or substitutes, if the ORIGINAL purchaser has completed & returned the Warranty Registration within 30 days of purchase & includes original receipt. Shipping, customs, duties & taxes must be PRE-PAID TO & FROM ObusForme by the PURCHASER. This warranty gives you rights that vary by province/state. This warranty may change.

#### WHAT IS NOT COVERED

Wear & tear, aging (including foam discoloration, flattening, density, consistency), accidental damages, alterations, mishandling, faulty adjustment, misuse, improper care, power damage, rental use, discontinued items, service by anyone other than ObusForme, shipping damages, over inflation, neglect, items sold 'as is' or damage due to natural acts are NOT covered.

**HOW TO OBTAIN WARRANTY SERVICE**  
Please obtain a Return Authorization Number & instructions prior to sending your item or it may be denied. Please inform Customer Service by:

Mail:  
HoMedics Group Canada  
6460 Kennedy Road, Unit C  
Mississauga, Ontario, Canada L5T 2X4  
Tel: (416) 785-1386 Fax: (416) 785-5862  
Toll Free: 1-888-225-7378  
8:30 a.m. to 5:00 p.m., Mon - Fri EST  
[www.obusforme.com](http://www.obusforme.com)

# OBUSFORME®

Corporate Member of ACE  
**ACE**

ASSOCIATION OF CANADIAN ERGONOMISTS

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