**OBUSFORME®**  
ENGINEERED FOR THE HUMAN BODY™

## SIDE TO SIDE WITH MASSAGE

SS-BLK-01

### Congratulations

on your purchase of this ObusForme Side to Side with Massage lumbar support cushion. This cushion gently supports the natural curvature of your lower back, preventing the tendency to slouch and helping to enhance overall sitting position, while the two speed massage function helps to sooth and relax your back muscles.

Notice a slight odor? Please note that some minor odor may exist with this product when it is first removed from the packaging. If so, we recommend that you allow the product to air out for a day or two prior to use.

### IMPORTANT SAFEGUARDS

**CAUTION: PLEASE READ ALL INSTRUCTIONS CAREFULLY BEFORE OPERATING.**

- This product is not intended for medical use. It is intended only to provide daily support and/or a soothing massage.
- Do not place or store this cushion where it can fall or be dropped into a tub or sink.
- Do not place or drop into water or any other liquid.
- Never operate this appliance if it has been damaged or if it is not working properly. Return it to the ObusForme for examination and repair.
- Keep away from heated surfaces.
- Do not use the massage function for more than 20 minutes at a time.

### HOW TO USE YOUR LUMBAR SUPPORT WITH MASSAGE

Small, lightweight and portable, your Lumbar Support with Massage is ideal for use anywhere you sit – at home, in the office and on the go.

### INSTALL BATTERIES

1. Unzip cover to expose battery compartment located in the cushion's foam. The zipper is located on the back of the cushion and is not designed to unzip completely, only to provide access to the battery compartment. (The cushion cover is not designed to be removed from the cushion.)
2. Gently pull the battery compartment away from the foam so that all four battery locations are exposed (two on the top, two on the bottom of the compartment).
3. Insert four new AA (LR6) alkaline batteries (not included). Take care to insert batteries as indicated inside the compartment.
4. Replace the battery compartment in the foam and close the zipper on the cover.

### BATTERY PRECAUTIONS

- Use only the size and type of batteries specified.
- When installing batteries, observe proper +/- polarities. Incorrect installation of batteries may cause damage to the unit.
- Do not mix old and new batteries.
- Do not mix alkaline, standard (carbon-zinc) or rechargeable (ni-cad, ni-mh, etc.) batteries.
- Do not dispose of batteries in fire. Batteries may explode or leak.
- If the unit is not going to be used for a long period of time, remove batteries to prevent damage due to possible battery leakage.
- Clean the battery contacts and also those of the device prior to battery installation.

- Remove depleted batteries promptly.
- Dispose of properly in accordance with all state, province and country regulations. This product contains small parts that may present a choking hazard for young children

### SITUATE CUSHION AND OPERATE THE MASSAGE FUNCTION

1. Simply place the cushion against the back of your chair or seat in an upright position, wherever it feels the most comfortable for you.
2. Use the attached adjustable elastic strap to secure it in place.
3. To turn massage on, adjust the switch located on the side of the cushion to "LO" for low intensity or "HI" for high intensity. To turn off, adjust the switch to the "OFF" position.

### USE & CARE TIPS

If the surface of the Side to Side cushion becomes soiled, it can be wiped gently with a damp cloth, mild soap and a sponge, or a soft brush. Do not rub the product excessively or machine wash. Do not soak.

### OBUSFORME® LIMITED WARRANTY

#### (CANADA Only)

ObusForme guarantees this product will be free from defects in workmanship & materials for a period of one (1) year from the date of original purchase. This guarantee applies when the product is used for the purpose intended. In the advent of a defect as described above, please return this product to ObusForme for refund or replacement. This warranty gives you rights that vary by province/state. This warranty may change.

#### WHAT IS NOT COVERED

Wear & tear, aging (including foam discoloration, flattening, density, consistency), accidental damages, alterations, mishandling, faulty adjustment, misuse, improper care, power damage, rental use, discontinued items, service by anyone other than ObusForme, shipping damages, neglect, items sold 'as is' or damage due to natural acts are NOT covered.

#### HOW TO OBTAIN WARRANTY SERVICE

Please obtain a Return Authorization Number & instructions prior to returning your item or it may be denied. Please contact Customer Service by:

Telephone: 1-888-225-7378 or 416-785-1386  
Mon – Fri, 8:30 a.m. to 5:00 p.m., EST

Fax: (416) 785-5862

Email: CustomerServiceCAN@obusforme.com

Mail:  
HoMedics Group Canada  
6460 Kennedy Road, Unit C  
Mississauga, Ontario, Canada L5T 2X4

JUN. 2017

Register your product (s) to ensure faster and more efficient customer service. It will allow us to keep you up-to-date with warranty and product information. Please take a few minutes to register online now.

Register at  
<http://www.obusforme.com/product-registration/>

**OBUSFORME®**

Corporate Member  
Membre Corporatif  
**ACE**

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Canadian  
Chiropractic  
Association

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