The ObusForme® Soft Gel Hot & Cold Compress XL

This reusable compress is designed for use with common aches, pains, swelling and bruising. Please read this carefully before use.

Cold Use

Recommended directly after an injury.

For optimum cold effectiveness, this compress should be stored in the freezer for a minimum of 2 hours. This compress will not freeze solid; it remains pliable to conform to the shape of the body.

When applying the compress, wrap it in a thin cloth that will not impede cold transfer. Cold therapy can be applied for 5-20 minutes at a time, depending on the location and size of the injury.

Hot Use

Not recommended for at least 2 days after an injury.

Stovetop: Immerse the compress in hot or boiling water for 5 to 7 minutes.

Microwave: Fully immerse the compress in water in a microwave-safe container. Heat on high power for 3 minutes. If further heating is required, microwave at 30-second intervals.

Safely remove the heated compress from the water and wrap it in a thin cloth that will not impede heat transfer. Heat can be applied for as long as 20-30 minutes. Comfort level and heat tolerance should be used as guidelines.

CAUTION!

If you have concerns regarding your health, please consult your physician before use. An injury that is really painful or pain that does not resolve itself after 24-48 hours of self-care should be checked by a professional. Do not apply the compress directly to the skin. To avoid burns, use the hot compress carefully. Use only as directed. Individuals with circulatory or diabetic problems should consult a physician before using this product. Do not ingest. Do not puncture. Contents may stain. Non-toxic – Safe for use with children.

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