## The Obus Forme® Soft Gel Hot & Cold Compress

This reusable compress is designed for use with common aches, pains, swelling and bruising. Please read this carefully before use.

### Cold Use

Recommended directly after an injury.

For optimum cold effectiveness, this compress should be stored in the freezer for a minimum of 2 hours. This compress will not freeze solid; it remains pliable to conform to the shape of the body.

When applying the compress, wrap it in a thin cloth that will not impede cold transfer. Cold therapy can be applied for 5-20 minutes at a time, depending on the location and size of the injury.

### Hot Use

Not recommended for at least 2 days after an injury.

Stovetop: Immerse the compress in hot or boiling water for 5 to 7 minutes.

Microwave: Fully immerse the compress in water in a microwave-safe container. Heat on high power for 3 minutes. If further heating is required, microwave at 30-second intervals.

Safely remove the heated compress from the water and wrap it in a thin cloth that will not impede heat transfer. Heat can be applied for as long as 20-30 minutes. Comfort level and heat tolerance should be used as guidelines.

### CAUTION!

If you have concerns regarding your health, please consult your physician before use. An injury that is really painful or pain that does not resolve itself after 24-48 hours of self-care should be checked by a professional. Do not apply the compress directly to the skin. To avoid burns, use the hot compress carefully. Use only as directed. Individuals with circulatory or diabetic problems should consult a physician before using this product. Do not ingest. Do not puncture. Contents may stain. Non-toxic – Safe for use with children.

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