# OBUSFORME®

## DOSSIER ULTRAFORME™

UF-SDM-LG, UF-SDM-MD, UF-SDM-SM,  
UF-BLK-SM, UF-BLK-MD, UF-BLK-LG, UF-OUS-CB

### CONGRATULATIONS

on your purchase of the UltraForme™ Backrest. Our unique convex and concave design helps align your spine into its proper curvature, as well as aiding in the relief and prevention of back problems.

### SCIENTIFICALLY TESTED

The UltraForme Backrest relieves back pain by offering complete spinal support while minimizing damaging pressure points and disc pressure. Properly used, it will also significantly reduce muscle strain in the lower back and increase respiratory capacity.

For optimal results, use the UltraForme Backrest with the ObusForme Contoured Seat Cushion to create a complete, comfortable and ergonomic sitting experience for your home, office or car. Together, these unique products encourage proper alignment between the pelvis and spine and further reduce pressure in the back, providing even more comfort for you each and every day.

### USING YOUR OBUSULTRAFORME BACKREST

1. Unfasten the buckles at the back of the backrest and place the backrest against the back of the seat/ chair where the device will be used.
2. Sit back and align your back against the UltraForme Backrest. Move the backrest up and down until the curvature of your UltraForme Backrest most comfortably meets your lumbar curvature (the small of your back).
3. Use the two quick-connect straps at the rear of the backrest to tightly secure the backrest to the chair in your ideal position.

4. Adjust the position of the back of your chair by tilting it backwards or forward as necessary.
5. Before using your UltraForme Backrest in a vehicle, please consult the owner's manual to ensure that you can properly adjust your vehicle's seat and/or headrest after installation of the UltraForme Universal Backrest.

Please Note: The UltraForme Backrest may in some cases require time to get used to. Over the years, your spine has become less flexible and needs time to adjust to its natural shape again. If you feel discomfort when using your UltraForme Backrest for the first time, it is recommended that you use it for shorter periods of time, gradually increasing your usage until you feel comfortable. If you have a very sore back or small curvature (lumbar area), it is highly recommended that you allow yourself a break-in period of at least 30 days in order to obtain the full benefits of your UltraForme Backrest. Those with scoliosis or a history of disc fusion surgery may not get the full benefit of the UltraForme Backrest. For these individuals, we recommend professional advice to determine if the ObusForme® Lowback Backrest Support is right for you.

### EASY CARE

Most cleaning can be easily accomplished with the use of a mild soap and water solution. However, we do suggest that you do a spot cleaning trial on an inconspicuous area. Do not rub the cover or immerse the entire cover in water as it may shrink.

### OBUSFORME LIMITED WARRANTY

ObusForme guarantees that the UltraForme will be free from defects in workmanship & materials for a period of one (1) year from the date of original purchase. This guarantee applies when the UltraForme is used for the purpose intended. In the event of a defect as described above, the UltraForme will be replaced with a new or refurbished product &/or substitute if the purchaser can provide the original purchase receipt. Shipping, customs, duties & taxes must be PRE-PAID TO & FROM ObusForme by the PURCHASER. This warranty gives you rights that vary by province/state. This warranty may change.

#### WHAT IS NOT COVERED

Wear & tear, aging (including foam discoloration, flattening, density, consistency), accidental damages, alterations, mishandling, faulty adjustment, misuse, improper care, power damage, rental use, discontinued items, service by anyone other than ObusForme, shipping damages, neglect, items sold 'as is' or damage due to natural acts are NOT covered.

#### HOW TO OBTAIN WARRANTY SERVICE

Please obtain a Return Authorization Number & instructions prior to sending your item or it may be denied. Please inform Customer Service by

Telephone: (416) 785-1386  
Toll Free: 1-888-225-7378  
8:30 a.m. to 5:00 p.m., Mon - Fri EST

HoMedics Group Canada  
6460 Kennedy Road, Unit C  
Mississauga, Ontario, Canada L5T 2X4

[www.ObusForme.com](http://www.ObusForme.com)

## OBUSFORME®

Corporate Member  
Membre Corporatif  
**ACE**  
Association of Canadian Ergonomists  
Association Canadienne d'Ergonomie

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HoMedics Group Canada  
6460 Kennedy Road. Bureau C  
Mississauga, Ontario, Canada L5T 2X4

[www.ObusForme.com](http://www.ObusForme.com)

**OBUSFORME**

Corporate Member  
Membre Corporatif  
**ACE**

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