PROFESSIONAL BODY MASSAGER™

Thank you for purchasing the Obus Forme Professional Body Massager. When using this product, there are basic precautions that should always be followed. Please read this notice carefully before use.

FEATURES

- 3" arm extension for hard-to-reach areas
- Variable speed control provides a gentle to exhilarating massage with the slide of a switch
- On/off switch
- Neoprene covered handle for a comfortable and firm grip
- 9.5' line power cord
- Three custom sets of massage head attachments offer soft or firm massage applications
- Dual pulsating heads deliver a kneading massage of up to 3,000 pulses per minute
- Convenient front handgrip helps to stabilize the unit for a concentrated massage
- Cord wrap strap with hook-and-loop enclosure (not shown)

IMPORTANT SAFETY INSTRUCTIONS

Warning! To reduce the risk of electric shock, fire, burns and/or injury:

- Use this product for its intended use only, as described in this manual. Do not use any attachments not recommended by Obus Forme.
- This product is designed for use in a 110V electrical outlet, in Canada and the U.S. only. Using this product with an electrical outlet of higher or foreign voltage will damage the unit and may cause fire, shock and/or injury. This appliance has a polarized plug (one blade is wider than the other).
- Do not operate unit under a blanket or pillow. Excessive heating can occur and cause fire, electric shock or injury.
- Always unplug the unit from the electrical outlet immediately after using, before cleaning and before attaching or removing parts.
- Do not use, place or store the unit where it can fall or be pulled into a tub or sink. Do not allow the unit to come into contact with water or any other liquid.
- Do not reach for the unit if it has fallen into water. Unplug it immediately. Do not operate.
- Do not open the unit to check for damages.
- Do not operate the unit if it has a damaged cord or plug, if it is not working properly or if it has been dropped or damaged in any manner. Contact Obus Forme for instructions for returning it for examination.
- Do not carry this appliance by its supply cord or use the cord as a handle.
- Keep cord away from heated surfaces.
- Ensure massage head attachments are secured onto unit before use.
- Do not block air openings or place unit on a soft surface, such as a bed or couch, where the air openings may be blocked. Keep air openings free of lint, hair, etc.
- Do not drop or insert foreign objects into any opening.
- Do not operate where aerosol (spray) products are being used or where oxygen is being administered.
- Do not operate the unit while in a sleepy state.
- Do not cover the unit while in operation.
- Do not use the unit directly on one concentrated area of the body. Vibrating massagers are designed to be used in a continuous motion, constantly changing the massage area.
- Do not use on or near eyes or other highly sensitive areas such as the head/skull, front of the neck, over the throat or chest, over or near the ears or groin.
- Keep small body parts such as fingers and hair away from massage heads and the gap between massage heads and the housing of the unit. This product was designed for use on broad body areas only.
- This appliance is designed for household use only; do not use outdoors.
- Close supervision is necessary when this appliance is used by, on or near children and disabled persons.

CAUTION!

- This product is not intended to replace the service or therapy provided by a professional or physician. Do not use this product as a substitute for medical attention.
- If you have any concerns regarding your health, consult your doctor before using this product.
- Do not use this product directly on swollen or inflamed areas, open wounds or skin eruptions.
- Do not use this product over the site of a tumor or mass.
- Do not use this product on limbs with impaired circulation (e.g. varicose veins, phlebitis, peripheral vascular disease).
- Do not use this product on the back on the neck if you have poor circulation to the brain (i.e. cerebral insufficiency, stroke).
- Individuals with epilepsy should not use this product.
- Individuals with osteoporosis should not use this product.
- Individuals who are taking blood thinners (anticoagulants), have a bleeding disorder or who bruise easily should not use this product.
- If you have a pacemaker, please consult a physician before using this product.
- If you are pregnant, please consult a physician before using this product.
- Do not use this product for more than 20 minutes at a time. Extensive use could lead to excessive heating and a shorter product life. Should the massager feel hot, discontinue used and allow the unit to cool before operating again.
- Do not use this product before bed, as the massage has a stimulating effect and can delay sleep.
- This product should not be used by individuals suffering from any physical ailment that would limit their capacity to operate the controls.

OPERATING YOUR PROFESSIONAL BODY MASSAGER

- Make sure the unit is set to the OFF position.
- Select a set of massage head attachments. Screw each attachment on clockwise. To change to a different set, simply unscrew, counter-clockwise, the set that is on the unit. Screw on the desired set.
Helpful hint: To further customize your massage, try using a different attachment on each head of the unit.
- Plug the unit into a 110-volt household electrical outlet (for use in Canada and U.S. only). This plug will fit in a polarized outlet only one way. If the plug does not fit fully in the outlet, reverse the plug. If it still does not fit, contact a qualified electrician to install a proper outlet. Do not change the plug in any way.
- Turn the unit on by sliding the ON/OFF switch to the ON position.
- Adjust the intensity of the massage by sliding the variable speed control switch upwards for a higher intensity and downwards for a lower intensity.
- Lightly apply the massage heads to the area you wish to treat using a slow, continuous motion.
- Treat area for a short period of time and then move to another area. Never treat the same area for longer than 3 minutes. Operation of this product for more than 20 minutes is not recommended and may cause overheating.
- To reach hard-to-reach areas, extend the arm of the unit by rotating the ring, located at the top of the arm, to the left. Gently pull the end of the handle and customize the extension at 1.5" or 3". To secure the length of the arm, simply rotate the ring to the right.
- Use the front handgrip for a concentrated, steady massage.
- When treatment is complete, turn the unit off by sliding the ON/OFF switch to the OFF position. Unplug the power cord and then wrap it up using the cord wrap strap.

CARING FOR YOUR PROFESSIONAL BODY MASSAGER

Servicing of the Professional Body Massager should only be performed by authorized Obus Forme personnel.

To Clean

- Unplug the unit and allow it cool before cleansing.
- Wipe using a soft, slightly damp cloth.
- Never allow water or any liquids to come into contact with the handle.

To Store

- Unplug the unit from the outlet and allow it to cool before storing it in a box or a clean, dry place.
- Wrap the power cord using the cord wrap strap. Do not wrap the power cord around the unit as this can cause the cord to wear and break.
- Do not hang the unit by its power cord.

LIMITED WARRANTY:

(North America Only) Obus Forme guarantees all items are free from defects in workmanship & materials for a time stated below from the original purchase date. This applies when items are used for the purpose intended. Items will be repaired/replaced with new/refurbished parts/items &/or alternates (our option) if the ORIGINAL purchaser has obtained Return Authorization (RA) from Customer Service and has sent the item along with its ORIGINAL receipt. Shipping, customs, duties and taxes must be PRE-PAID TO and FROM OBUS FORME by the PURCHASER.

This warranty gives you rights that vary by province/state. This warranty may change. This warranty is non-transferable.

WHAT IS NOT COVERED

Wear & tear, aging, foam/item discoloring, odor, flattening, density, variation, leaking, alteration, mishandling, faulty adjustment, misuse, improper care, power damage, accidents, rental use, obsolete items, service by anyone other than Obus Forme, use of any non-Obus Forme authorized parts, shipping damage, over-inflation, neglect, items sold 'as is' or damage due to natural acts are NOT covered.

WARRANTY TIME FRAME

Backrest Frame: Lifetime (cover, foam, lumbar pad and other Backrest parts/materials are NOT covered)
Seat Frame: 1 year (cover, foam and other Seat parts/materials are NOT covered)
Back Therapy: 1 year (includes but is NOT limited to Back Belts, Backpacks, Drivers Seats)
Sleep/Foot/Muscle Therapy: 1 year (pillow case is NOT included & must not be returned. Only foam or bladder will be replaced.)
Electrical Parts: 1 year (This includes wires, adaptors, plugs and other electrical parts/components)

HOW TO OBTAIN WARRANTY SERVICE

You must obtain Return Authorization & direction before sending your item or it will be DENIED. Please obtain an RA via Customer Service by:

Mail:

HoMedics Group Canada
 344 Consumers Road, Toronto, Ontario,
 Canada M2J 1P8
 Tel: (416) 785-1386 Fax: (416) 785-5862
 Toll Free: 1-888-225-7378
 8:30 a.m. to 5:00 p.m., Mon - Fri ET
 www.obusforme.com

ASSOCIATION OF CANADIAN ERGONOMISTS

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Obus Forme is committed to providing you with optimal relief and comfort. To serve you better in the future, we would like to know if we have fulfilled our commitment. Please complete and return this Questionnaire to help us better meet your needs.

We aggregate this information and use it internally for research and marketing purposes only. We do not disclose personal information to any third parties. If you have any questions about the personal information that we keep on file, please contact a customer service representative at the number listed below.

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Age / Âge : _____

Occupation / Profession : _____

1. Which Obus Forme® product did you purchase? / Quel produit Obus Forme® avez-vous acheté?

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